## **SPRINT - 3**

Date	12 November 2022		
Team ID	PNT2022TMID42447		
Project Name	Personal Assistance for Seniors Who Are Self- Reliant.		

## 1. Python code for random medicine and time generating:

```
import json
import wiotp.sdk.device
import time
import random
myConfig = {
  "identity": {
    "orgId": "mni3qc",
    "typeId": "medicine",
    "deviceId": "123456"
   "auth": {
    "token": "paul@123"
client = wiotp.sdk.device.DeviceClient(config=myConfig, logHandlers=None)
client.connect()
for i in range(0,20):
  tablet=["Paracetamol","Aspirine","Azithral","Asthalin","Sinarest"]
  medicinetime=[12.00,1.00,2.00,3.00,5.00,18.00,20.00,7.00]
  name = "mani"
  medicine=random.choice(tablet)
  medicinetime=random.choice(medicinetime)
  mydata = {'Patient Name': name, 'Medicine Name': medicine, 'Time': medicinetime}
  client.publishEvent("IoTSensor", "json", data=mydata, qos=0, onPublish=None)
  print("Data published to IBM IOT platform :", mydata)
  time.sleep(5)
client.disconnect()
```

```
temp.py - C:\Users\91637\Downloads\temp.py (3.7.4)
```

```
- 0 ×
File Edit Format Run Options Window Help
```

```
wiotp.sdk.device
            time
import random
myConfig = {
       "identity": {
    "orgId": "mni3qc",
    "typeId": "medicine",
    "deviceId": "123456"
          "auth": {
    "token": "paul@123"
      }
client = wiotp.sdk.device.DeviceClient(config=myConfig, logHandlers=None)
medicine=random.choice(tablet)
| medicinerrandom.choice(tablet)
medicinerimerrandom.choice(medicinetime)
medicinetime:random.choice(medicinetime)
mydata = {'Fatient Name': name, 'Medicine Name': medicine, 'Time': medicinetime)
client.publishEvent("foTSensor", "json", data-mydata, qos=0, onFublish=None)
print("Data published to IBM IOT platform :", mydata)
time.sleep(5)
client.disconnect()
```

## \*Python 3.7.4 Shell\* - n ×

```
File Edit Shell Debug Options Window Help
Python 3.7.4 (tags/v3.7.4:e09359112e, Jul 8 2019, 20:34:20) [MSC v.1916 64 bit (AMD64)] on win32
Type "help", "copyright", "credits" or "license()" for more information.
2022-11-16 15:50:06,267 wiotp.sdk.device.client.DeviceClient INFO Connected successfully: d:mni3qc:medicine:123456
Data published to IBM IOT platform: ('Patient Name': 'mani', 'Medicine Name': 'Asthalin', 'Time': 5.0)
Data published to IBM IOT platform: ('Patient Name': 'mani', 'Medicine Name': 'Asthalin', 'Time': 12.0)
Data published to IBM IOT platform: ('Patient Name': 'mani', 'Medicine Name': 'Asithral', 'Time': 12.0)
Data published to IBM IOT platform: ('Patient Name': 'mani', 'Medicine Name': 'Asithral', 'Time': 18.0)
Data published to IBM IOT platform: ('Patient Name': 'mani', 'Medicine Name': 'Aspirine', 'Time': 18.0)
Data published to IBM IOT platform: ('Patient Name': 'mani', 'Medicine Name': 'Aspirine', 'Time': 18.0)
Data published to IBM IOT platform: ('Patient Name': 'mani', 'Medicine Name': 'Aspirine', 'Time': 20.0)
Data published to IBM IOT platform: ('Patient Name': 'mani', 'Medicine Name': 'Asithral', 'Time': 12.0)
Data published to IBM IOT platform: ('Patient Name': 'mani', 'Medicine Name': 'Asithral', 'Time': 12.0)
Data published to IBM IOT platform: ('Patient Name': 'mani', 'Medicine Name': 'Sinarest', 'Time': 1.0)
Data published to IBM IOT platform: ('Patient Name': 'mani', 'Medicine Name': 'Asthalin', 'Time': 12.0)
Data published to IBM IOT platform: ('Patient Name': 'mani', 'Medicine Name': 'Asthalin', 'Time': 5.0)
Data published to IBM IOT platform: ('Patient Name': 'mani', 'Medicine Name': 'Sinarest', 'Time': 1.0)
Data published to IBM IOT platform: ('Patient Name': 'mani', 'Medicine Name': 'Sinarest', 'Time': 1.0)
Data published to IBM IOT platform: ('Patient Name': 'mani', 'Medicine Name': 'Sinarest', 'Time': 1.0)
Data published to IBM IOT platform: ('Patient Name': 'mani', 'Medicine Name': 'Sinarest', 'Time': 2.0)
Data published to IBM IOT platform: ('Patient Name': 'mani', 'Medicine Name': 'Sinarest', 'Time': 2.0)
Data published to IBM IOT platform: ('Patient Name': 'mani', 'Medicine Name': 'Sinarest', 'Time': 2.0)
Data published to IBM IOT platform: ('Patient Name': 'mani', 'Medici
                                                 ----- RESTART: C:\Users\91637\Downloads\temp.pv ----
```

